

# Solvent systems - Lipid migration

Lipids	Solvent system (Rf) [see chart for system]					
	i	ii	iii	iv	v	vi
Triglycerides					0.70	0.96
Free Fatty Acid					0.51	0.16
Diglycerides 1,2					0.70	0.24
Diglycerides 1,3						0.32
Monoglycerides						
Phosphatidylethanolamine	0.79	0.55		0.43		
Phosphatidyl (Monomethylethanolamine)	0.71	0.41		0.33		
Cardiolipin	0.67	0.56	0.38			
Phosphatidylglycerol	0.60	0.50	0.31	0.36		
Phosphatidyl (Dimethylethanolamine)	0.58	0.56		0.27		
Phosphatidic Acid	0.55	0.05	0.58			
Phosphatidylinositol	0.39	0.10				
Phosphatidylcholine	0.34	0.30		0.17		
Phosphatidylserine	0.33	0.12				
Cerebrosides	0.94	0.55				
Ceramides	sf*	sf*				
Sphingosine	0.28	0.75				
Sphingomyelin	0.28	0.13				
Synthetic Cardiolipin	0.29	0.20				
Lyso-Phosphatidylglycerol	0.54	0.20				
Lyso-Phosphatidylethanolamine	0.45	0.20				
Monolysocardiolipin	0.45	0.34				
Lyso-Phosphatidic Acid	0.40	0.01				
Dilysocardiolipin	0.32	0.21				
Lyso-Phosphatidylinositol	0.29	0.03				
Lyso-Phosphatidylcholine	0.22	0.08				
Lyso-Phosphatidylserine	0.18	0.02				

Solvent systems		Ratio (V:V)
i	Chloroform:Methanol:Water	65:25:4
ii	Chloroform:Methanol:Ammonium Hydroxide	65:25:4
iii	Chloroform:Hexane:Methanol:Acetic Acid	50:30:10:5
iv	Toluene:Pyridine:Water	60:60:10
v	Cyclohexane:Ethyl Acetate	3:2
vi	Toluene:Chloroform:Methanol	85:15:5

### Notes

R<sub>f</sub> values represent relative migration only, whereas absolute values depend on various environmental parameters (e.g. temperature, humidity) which may vary depending on location.

\*sf solvent front

### Non-warranty

The information in this publication is believed to be accurate and is given in good faith, but no representation or warranty as to its completeness or accuracy is made. Suggestions for uses or applications are only opinions. Users are responsible for determining the suitability of these products for their own particular purpose. No representation or warranty, expressed or implied, is made with respect to information or products including, without limitation, warranties of merchantability, fitness for a particular purpose, non-infringement of any third-party patent or other intellectual property rights including, without limit, copyright, trademark and designs. Unless otherwise stated, any trademarks identified herein are trademarks of the Croda group of companies. Avanti Research is a trademark of Croda international Plc. Avanti Research is a Croda brand associated with Avanti Polar Lipids, LLC.