



Lipids for liposome formation

Name	Fatty acid	Transition temperature (°C)	Net charge at pH 7.4	Product number
DLPC	12:0	-1	0	850335
DMPC	14:0	23	0	850345
DPPC	16:0	41	0	850355
DSPC	18:0	55	0	850365
DOPC	18:1	-20	0	850375
DMPE	14:0	50	0	850745
DPPE	16:0	63	0	850705
DOPE	18:1	-16	0	850725
DMPA•Na	14:0	50	-1.3	830845
DPPA•Na	16:0	67	-1.3	830855
DOPA•Na	18:1	-8	-1.3	840875
DMPG•Na	14:0	23	-1	840445
DPPG•Na	16:0	41	-1	840445
DOPG•Na	18:1	-18	-1	840475
DMPS•Na	14:0	35	-1	840033
DPPS•Na	16:0	54	-1	840037
DOPS•Na	18:1	-11	-1	840035
DOPE-Glutaryl•(Na) ₂	18:1	~-10	-2	870242
Tetramyristoyl Cardiolipin•(Na) ₂	14:0	59	-2	710332
DSPE-mPEG-2000•Na	18:0	N/A	-1	880120
DSPE-mPEG-5000•Na	18:0	N/A	-1	880220
DSPE-Maleimide PEG2000•Na	18:0	N/A	-1	880126
DOTAP•Cl	18:1	~-0	+1	890890

References

Marsh, D. (2012) Thermodynamics of phospholipid self-assembly. Biophys J 102:1079-87. [PubMed]

Non-warranty

The information in this publication is believed to be accurate and is given in good faith, but no representation or warranty as to its completeness or accuracy is made. Suggestions for uses or applications are only opinions. Users are responsible for determining the suitability of these products for their own particular purpose. No representation or warranty, expressed or implied, is made with respect to information or products including, without limitation, warranties of merchantability, fitness for a particular purpose, non-infringement of any third-party patent or other intellectual property rights including, without limit, copyright, trademark and designs. Unless otherwise stated, any trademarks identified herein are trademarks of the Croda group of companies. Avanti Research is a trademark of Croda international Plc. Avanti Research is a Croda brand associated with Avanti Polar Lipids, LLC.